



BioMEDICAL
LIFE-SYSTEMS

KIVUNLIEVITYKSEEN

IMPULSE TENS D5

- Digitaalinen 2- kanavainen tns
- 4 valmista ohjelmaa
- Yksilöllinen säätömahdollisuus 1-120Hz ja impulssin pituus 25-250μs
- Jatkuva virta, moduloitu taajuus, moduloitu pulssin pituus, burst I ja burst II
- Ajan säätö/jatkuva aika
- Toimii 2:lla sormipatterilla

Tekniset tiedot:

Size: 4.2" x 2.9" x 1.4" (10.7cm x 7.4cm x 4cm)

Weight: 3.35 oz (94.97 grams)

Power Source: 2 AA batteries, Type LR6

Channels: Dual

Waveform: Symmetrical, biphasic square wave

Pulse Rate (Hz): 1-120 Hz (Hertz or pps) adj.

Pulse Width (μS): 25-250 microseconds (μS) adjustable

Output: Constant current

Constant: Continuous stimulation. Pulse rate/Pulse Width are adjustable

Pulse Rate Modulation:

Pulse Rate modulates from 100Hz down to 20Hz over a 15 second cycle (7.5 seconds down, 7.5 seconds up). Pulse Width is adjustable.

Pulse Width Modulation:

Pulse width modulates from 125-250μS and back down over a 5 second cycle (2.5 seconds down, 2.5 seconds up). Pulse Rate is adjustable.

Burst I: One second ON. One second OFF. Pulse Rate and Pulse Width are adjustable.

Burst II: Seven Pulses per burst, 2 bursts per second. Pulse Width is adjustable Output: Constant current

Intensity: Continuously adjustable from 0-100 mA peak to peak

Output Voltage: Continuously adjustable from 0-50 V peak to peak

Tolerances: +/-1%